

1M10 Adjust

USA 800.328.4058
Fax 800.962.2549

Order Form

OB Customer No. _____

City/State/Zip _____

PO No. _____

Ordered By _____

Company _____

Phone Number _____ Date _____

Address _____

Shipping Options:
 Ground 2-Day Next Day p.m. Next Day a.m.
 Other _____

Email _____

1M10 Adjust

Quantity	Order no.	Side	Size	Weight category	Color	Shape
	1M10=			-	-P/	
	1M10=			-	-P/	

Order example
1M10=L27-3-P/4N

	Side
Right	R
Left	L

	Size (cm)
Normal (N)	22, 23, 24, 25, 26, 27, 28, 29, 30
Small (S)	22, 23, 24, 25, 26

	Color
Beige	4
Brown	15

	Foot shape
Normal foot shape, heel height 10 mm (3/8")	N
Small foot shape, heel height 20 mm (3/4")	S

Size (cm)				Weight category
22-23	24-25	26-27	28-30	
to 52 kg (to 115 lbs)	to 58 kg (to 128 lbs)	to 72 kg (to 158 lbs)	to 77 kg (to 170 lbs)	1
53 – 68 kg (116 – 150 lbs)	59 – 76 kg (129 – 169 lbs)	73 – 95 kg (161 – 210 lbs)	78 – 100 kg (172 – 220 lbs)	2
69 – 80 kg (151 – 176 lbs)	77 – 100 kg (169 – 220 lbs)	96 – 125 kg (211 – 275 lbs)	101 – 125 kg (221 – 275 lbs)	3



When ordering a 1M10 Adjust, the delivery will include the 2C1 Footshell.

Additional Components

Quantity	Article number

Quantity	Article number

Distribution Center: 14630 28th Avenue North, Minneapolis, MN 55447
 Manufacturing: 3820 West Great Lakes Drive, Salt Lake City, UT 84120
 For more information please visit our web site at www.ottobockus.com

All specifications and pricing are subject to change without notice.
 ©2010 Otto Bock HealthCare LP • 09091703.9 • 2/10

Otto Bock®

QUALITY FOR LIFE