

Otto Bock®

QUALITY FOR LIFE



KIWI Mobility System

Low Weight, High Design

For days when the to-do list is too long, the KIWI pediatric wheelchair makes daily travels smoother for parents and kids. At just 24 lbs, the KIWI helps make walking around the neighborhood and getting around the city a little easier.

In addition to its good looks, the KIWI helps connect parent and child with a reversible seat and handle bar. The resulting configurations let caregivers interact with the child while positioning the wheels to optimally cover different terrain. Additional adjustments come from the height-adjustable handle bar and seat.

Product Information

Features and Benefits

- Pivoting handlebar and reversible seat orientation allows child to have face-to-face interaction with the caregiver — or to see the world glide by
- 30° tilt and 30° recline offer pressure relief and positioning options to help calm the child
- Contemporary design adds style and mild positioning
- Spring loaded caster forks deliver a smoother ride
- Hand-operated wheel lock is easy to engage

Options

To further customize KIWI select from the following:

- › Five point harness
- › Removable grab rail
- › Abductor
- › Backrest positioning kit
- › Sun canopy
- › Storage bag
- › Available in two sizes



Central Unit



Grab Rail



Footplate



Pivoting the Handle Bar



Swivel Lock



Sun Canopy

Technical Data

	Size 1	Size 2
Seat Width	10.6" (27 cm)	12.6" (32 cm)
Seat Depth	4.7 – 8.7" (12 – 22 cm)	7.9 – 11" (20 – 28 cm)
Back Height	24.8" (63 cm)	26.8" (68 cm)
Handle Bar Height	38 – 48" (96 – 116 cm)	38 – 48" (96 – 116 cm)
Lower Leg Length	5.5 – 7.9" (14 – 20 cm)	7.7 – 11.6" (19.5 – 29.5 cm)
Tilt	0° – 30°	0° – 30°
Seat to Back Angle	0° – 30°	0° – 30°
Overall Width	23" (58 cm)	25.6" (65 cm)
Maximum Load Capacity	55 lbs (25 kg)	55 lbs (25 kg)
Weight (including seat unit)	24 lbs (11 kg)	26 lbs (12 kg)
Folding Size (including seat unit)	32" x 23" x 15" (80 x 58 x 37 cm)	32" x 25" x 20" (80 x 63 x 50 cm)

Otto Bock®

QUALITY FOR LIFE