

Date _____ Account Name _____ Bill to Address _____ City/State/Zip _____ Fax No. _____ Physician/Therapist _____	Account No. _____ Ship to Address _____ City/State/Zip _____ Phone No. _____ Email Address _____ Authorized Signature _____
---	--

Mark For: Quote Only Please Order

Purchase Order No.: _____

I. Nurmi NeoGait Trainer

For children requiring positioning aids and support during ambulation and various stages of gait development.

Includes:

- Reverse-roll locks on the rear wheels
- Swivel front wheels
- Puncture proof tires
- Foldable for easy transport
- 5 year limited warranty

Indicate Color:

- Saturn Blue Metallic – 306
 Sparkle Green – 311

			HCPCS
<input type="checkbox"/> Size 1	HR24160000	\$530.00	E8000
<input type="checkbox"/> Size 2	HR24260000	\$541.00	E8000
<input type="checkbox"/> Size 3	HR24360000	\$565.00	E8000


Please choose Grip Bar from below.

II. Postural Supports and Accessories

Standard Grip Bar

Continuous height and depth adjustable grip bar.


- Size 1 HR24161600 \$111.00
 Size 2 HR24261600 \$111.00
 Size 3 HR24361600 \$111.00



Grip Bar with Universal Grips

Grips can be adjusted to reduce the width between grips, as well as the grip height.

- Size 1 HR24161700 \$286.00
 Size 2 HR24261700 \$286.00
 Size 3 HR24361700 \$286.00




II. Postural Supports and Accessories (continued)

Extra High Grip Bar with Forearm Supports and Vertical Hand Grips

Height, depth and angle adjustable


- Size 1 HR24161800 \$426.00
 Size 2 HR24261800 \$426.00
 Size 3 HR24361800 \$426.00



Extra High Grip Bar with Standard Grips

Continuous height and depth adjustable Grip Bar that allows for extra heights.


- Size 1 HR24161880 \$223.00
 Size 2 HR24261880 \$223.00
 Size 3 HR24361880 \$223.00



Mini Grip Bar with Universal Grips

When a lower grip height and/or shorter width between grips is needed for smaller children.


- Size One Only HR24168600 \$221.00



Universal Grip (each)

Sold individually as needed for grip height discrepancies


- One Size HR24161730E \$71.00



Adjustable Forearm Support complete (each)

Sold individually as needed for grip strength discrepancies

- Size 1 HR24161850E \$164.00
 Size 2 & 3 HR24261850E \$164.00



II. Postural Supports and Accessories (continued)

Hand Plate for Grip Bars

(Available for Standard and Universal grip bars only)

Prevents the child's hands from slipping forward off the end of the grip

- For Standard Grip Bar HR24160070 \$156.00
- For Universal Grip Bar HR24160075 \$156.00

Gluteal Support

To assist with maintaining upright trunk position. Height, depth and angle adjustable

- One Size HR24168300 \$194.00

Combined Gluteal and Pelvic Support (3D Adjustable)

To assist with maintaining upright trunk position. Height, depth and angle adjustable

- One Size HR24168400 \$333.00

Anti-Tip Bars (pair) (indicate color)

Provides stability to prevent posterior tipping.

- Saturn Blue Metallic – 306
- Sparkle Green – 311
- Size 1 HR24167300-* \$72.00
- Size 2 & 3 HR24267300-* \$72.00

Friction Brake

Applies friction to the rear wheels to slow forward progression/speed during ambulation. Friction can be increased or decreased by adjusting the interior set screws.

- One Size HR24167200 \$39.00

Caster Swivel Wheel Locks*

**Must be installed by manufacturer at time of order*

Caster wheels can be locked in straight alignment to help control lateral movement during ambulation. Quick release feature allows for quick unlocking as gait development progresses.

- One Size HR24167100 \$111.00

*Indicate Color

II. Postural Supports and Accessories (continued)

Fold Up Seat (indicate color)

Helps build endurance levels for ambulation over greater distances.

- Saturn Blue Metallic – 306
- Sparkle Green - 311
- Size 1 HR24167500-* \$190.00
- Size 2 HR24267500-* \$190.00
- Size 3 HR24367500-* \$190.00

Storage Basket, black

Allows personal items to be kept within reach.

- One Size HR24168200 \$86.00

*Mounted on leftside of grip bar, please request alternate location if needed.

Sling Seat

Helps build endurance ambulation over greater distances. Supports the child's weight when needed.

- Size 1 HR24167599 \$117.00
- Size 2 HR24267599 \$117.00
- Size 3 HR24367599 \$117.00

Drive Break and Wheel Lock

Helps the child control the speed and secure the walking aid on incline surfaces.

For Standard Grip Bar

- Size 1 HR24168100 \$175.00
- Size 2 HR24268100 \$175.00
- Size 3 HR24368100 \$175.00

For Universal Grip Bar

- Size 1 HR24168800 \$180.00
- Size 2 HR24268800 \$180.00
- Size 3 HR24368800 \$180.00

For Extra High Grip Bar

- Size 1 HR24168700 \$180.00
- Size 2 HR24268700 \$180.00
- Size 3 HR24368700 \$180.00

II. Nurmi Neo Technical Specifications

Dimensions (in) and Weight (lbs)	Size 1	Size 2	Size 3
Minimum -Maximum Handle Height			
Standard Grip Bar	17.5" - 22.4"	22.3" - 28.5"	25" - 33.5"
With forearm supports (measured to lowest pad surface) vertical grips	20" - 24.5"	25.25" - 30.5"	27" - 33"
Mini Grip Bar	12" - 19"	N/A	N/A
Universal Grip Bar	15" - 22.4"	20" - 28.3"	22.4" - 33.5"
Extra High Grip Bar with Vertical Grips			
To Vertical Grip	22.5" - 28.3"	25.5" - 34"	27.5" - 36"
To Forearm Supports (measured to lowest pad surface)	22.5" - 28.4"	25.5" - 34"	27.3" - 35.8"
Extra High Grip Bar with Standard Grips	27"	32"	33.75"
Inside Width between Handlebars			
Distance between Grips: Standard Grip Bar	15"	17"	19"
Distance between Grips: Mini Grip Bar	8" - 15"		
Distance between Grips: Universal Grip Bar	10.6" - 15"	12.5" - 17"	14.5" - 19"
Distance between Forearm Supports (center)	9.8" - 15"	11.8" - 17"	13.8" - 19"
Forearm Support (width x length)(inside)	2.75" - 6.25"	2.75 - 8.25"	2.75 - 8.25"
Tube Diameter			
Grip	1"	1"	1"
Frame	1.25"	1.25"	1.25"
Overall Width	25"	26.7"	28"
Overall Length (with anti-tipper add 4")	25.5"	30"	32.5"
Folding Dimensions (L x H x W)			
Length	25.5"	29.5"	33.5"
Height	24.8"	26.8"	28.3"
Width	9.5"	9.5"	9.5"
Turning Diameter (without accessories)	34"	39"	42"
Seat Height	11.7"	17.6"	18.7"
Front Wheel Diameter	6"	6"	6"
Rear Wheel Diameter	6"	8"	8"
Weight Capacity	55 lbs.	88 lbs.	121 lbs.
Product Weight	12 lbs.	13 lbs.	14 lbs.

Since we are constantly updating and adding to our line of products, certain items may be discontinued or redesigned without notice.