

Retail Price List / Order Form

Price Effective January 18/2010 Rev. 4

Date	Ship To
Bill To	Address
Account #	City/Prov/PC
Address	Fax No.
Phone No.	Ordered By
E-Mail Address	Authorized Signature
Contact Person	

<input type="checkbox"/> Quote Only	<input type="checkbox"/> Please Order	Purchase Order No. _____
-------------------------------------	---------------------------------------	--------------------------

I. Mobility Base

For children requiring positioning aids and support during ambulation and various stages of gait development includes:

- Reverse-roll locks on the rear wheels
- Swivel front wheels
- Puncture proof tires
- Foldable for easy transport
- 5 year limited warranty

Indicate Colour:

- Saturn Blue Metallic – 306
 Granny Smith Green – 311

<input type="checkbox"/> Size 1	HR24160000	\$901.00
<input type="checkbox"/> Size 2	HR24260000	\$913.00
<input type="checkbox"/> Size 3	HR24360000	\$933.00

II. Postural Supports and Accessories

Standard Grips

<input type="checkbox"/> Size 1	HR24161600	\$189.00
<input type="checkbox"/> Size 2	HR24261600	\$189.00
<input type="checkbox"/> Size 3	HR24361600	\$189.00

Grip Bar with Universal Grips

Grips can be adjusted to reduce the width between grips, as well as the grip height.

<input type="checkbox"/> Size 1	HR24161700	\$248.00
<input type="checkbox"/> Size 2	HR24261700	\$272.00
<input type="checkbox"/> Size 3	HR24361700	\$309.00

Extra High Grip Bar with Forearm Supports + Vertical Hand Grips

Height, depth and angle adjustable.

<input type="checkbox"/> Size 1	HR24161800	\$435.00
<input type="checkbox"/> Size 2	HR24261800	\$435.00
<input type="checkbox"/> Size 3	HR24361800	\$435.00

II. Postural Supports and Accessories (continued)

Combined Gluteal and Pelvic Support (3D Adjustable)

To assist with maintaining upright trunk position. Height, depth and angle adjustable.

<input type="checkbox"/> One Size	HR24168400	\$576.00
-----------------------------------	------------	----------

Anti-Tip Bars (pair)

Provides stability to prevent posterior tipping.

Indicate Colour:

- Saturn Blue Metallic – 306
 Granny Smith Green – 311

<input type="checkbox"/> Size 1	HR24167300	\$163.00
<input type="checkbox"/> Size 2 & 3	HR24267300	\$163.00

Friction Brake

Applies friction to the rear wheels to slow forward progression/speed during ambulation. Friction can be increased or decreased by adjusting the interior set screws.

<input type="checkbox"/> One Size	HR24167200	\$121.00
-----------------------------------	------------	----------

Caster Swivel Wheel Locks*

**Must be installed by manufacturer at time of order*

Caster wheels can be locked in straight alignment to help control lateral movement during ambulation. Quick release feature allows for a quick unlocking as gait development progresses.

<input type="checkbox"/> One Size	HR24167100	\$204.00
-----------------------------------	------------	----------

Fold Up Seat

Helps build endurance levels for ambulation over greater distances.

Indicate Colour:

- Saturn Blue Metallic – 306
 Granny Smith Green – 311

<input type="checkbox"/> Size 1	HR24167500	\$335.00
<input type="checkbox"/> Size 2	HR24267500	\$335.00
<input type="checkbox"/> Size 3	HR24367500	\$335.00

*All specifications and pricing are subject to change without notice

Additional Options

Nurmi Neo

<input type="checkbox"/> Basket	HR24168200	\$171.00
<input type="checkbox"/> Sling Seat, Size 1	HR24167599	\$262.00
<input type="checkbox"/> Sling Seat, Size 2	HR24267599	\$262.00
<input type="checkbox"/> Sling Seat, Size 3	HR24367599	\$262.00
<input type="checkbox"/> Hand Grip	HR24160070	\$214.00
<input type="checkbox"/> Hand Grip 22cm	HR24160075	\$214.00
<input type="checkbox"/> Driving Brake for Standard Grips, Sz 1	HR24168100	\$205.00
<input type="checkbox"/> Driving Brake for Extra High Grips, Sz 1	HR24168700	\$248.00
<input type="checkbox"/> Driving Brake for Universal Grips, Sz 1	HR24168800	\$248.00
<input type="checkbox"/> Driving Brake for Standard Grips, Sz 2	HR24268100	\$205.00
<input type="checkbox"/> Driving Brake for Extra High Grips, Sz 2	HR24268700	\$248.00
<input type="checkbox"/> Driving Brake for Universal Grips, Sz 2	HR24268800	\$248.00
<input type="checkbox"/> Driving Brake for Standard Grips, Sz 3	HR24368100	\$205.00
<input type="checkbox"/> Driving Brake for Extra High Grips, Sz 3	HR24368700	\$248.00
<input type="checkbox"/> Driving Brake for Universal Grips, Sz3	HR24368800	\$248.00

III. Nurmi Neo Technical Specifications

Dimensions (in)(cm) and Weight (lb)(kg)	Size 1	Size 2	Size 3
Minimum Handle Height:			
Standard Grip Bar	18" (45.7 cm)	22.5" (57.2 cm)	25" (63.5 cm)
Mini Grip Bar	12" (30.5 cm)	N/A	N/A
Maximum Handle Height:			
Standard Grip Bar	22" (55.9 cm)	30" (76.2 cm)	33.5" (85.1 cm)
Inside Width between Handlebars:			
Standard Grip Bar	16" (40.6 cm)	18" (45.7 cm)	20" (50.3 cm)
Mini Grip Bar	8"-16" (20.3cm-40.6 cm)	N/A	N/A
Universal Grip Bar	10"-16" (26.7cm-40.6 cm)	12.5"-18" (31.8 cm-45.7 cm)	14.5"-20" (36.8 cm-50.8 cm)
Overall Length (with anti-tipper add 4")	25.5" (64.8 cm)	30" (76.2 cm)	32.5" (82.6 cm)
Overall Width	25" (63.5 cm)	26.75" (67.9 cm)	28" (71.1 cm)
Weight Capacity	55 lb. (24.9 kg)	88 lb. (39.9 kg)	121lb. (54.9 kg)
Product Weight	12 lb. (5.4 kg)	12 lb. (5.0 kg)	14lb. (6.4 kg)
Front Wheel Diameter	6" (15.2 cm)	6" (15.2 cm)	6" (15.2 cm)
Rear Wheel Diameter	6" (15.2 cm)	8" (20.4 cm)	8" (20.4 cm)
Folding Diameter			
Length	25.5" (64.8 cm)	29.5" (74.9 cm)	33.5" (85.1 cm)
Width	25" (63.5 cm)	26.75" (67.9 cm)	28.25" (71.8 cm)
Depth	9.5" (24.1 cm)	9.5" (24.1 cm)	9.5" (24.1 cm)
Turning diameter (without accessories)	34" (86.4 cm)	39" (99.1 cm)	42" (106.7 cm)
Seat Height	11.7" (29.7 cm)	17.6" (44.7 cm)	18.7" (47.5 cm)

Since we are constantly updating and adding to our line of products, certain items may be discontinued or redesigned without notice.