

Models

Four models with different connections are available.

Torsion Adapter with Clamp (Ø 30 mm) and Pyramid Receiver:

The 4R85 Torsion Adapter connects the knee joint with the Ø 30 mm tube adapter. It replaces the tube clamp adapter.

Torsion Adapter with Clamp (Ø 34 mm) and Pyramid Receiver:

The 4R86 Torsion Adapter connects the knee joint with the Ø 34 mm tube adapter. It replaces the tube clamp adapter.

Torsion adapter with Tube and Pyramid Receiver:

The 4R39 Torsion Adapter contains a Ø 30 mm tube. It replaces the tube adapter and bridges the distance between prosthetic foot and knee joint.

Torsion Adapter with Screw Plate and Pyramid Receiver:

The 4R40 Torsion Adapter has a screw plate at its proximal end. It replaces the socket attachment block.

Torsion Adapter


Stainless steel

Torsion Adapter with Pyramid Receiver and Ø 30 mm Tube Connection.



Order no.	4R85
Diameter	30 mm
Weight	350 g
System height	68 mm
Material	Stainless steel
Mobility Grade	1 - 4
Max. body weight	100 kg

- ▶ When using the 4R85 with the knee joints with friction brake 3R15 and 3R49, the longer extension assist pull (included with torsion adapter) must be mounted.

 647G23




Torsion Adapter

Titanium

Torsion Adapter with Pyramid Receiver and Ø 34 mm Tube Connection.



Order no.	4R86
Diameter	34 mm
Weight	340 g
System height	68 mm
Material	Titanium
Mobility Grade	1 - 4
Max. body weight	110 kg

 647G23



≤ 110 kg

Torsion Adapter

Stainless steel

Torsion Adapter with Proximal Screw Plate and Distal Pyramid Receiver.



Order no.	4R40
Weight	340 g
System height	58 mm
Material	Stainless steel
Mobility Grade	1 - 4
Max. body weight	125 kg

647G23

Torsion Adapter

Stainless steel

Torsion Adapter with Ø 30 mm Tube and Pyramid Receiver.



Order no.	4R39
Diameter	30 mm
Weight	500 g
Min. system height	113 mm
Max. system height	476 mm
Material	Stainless steel
Mobility Grade	1 - 4
Max. body weight	100 kg

647G23

Single components

replacement parts

Order no.	4R39	4R40	4R85	4R86
506G3=M8x12 Set Screw	▲	▲	▲	▲
4D4 Single component Pack			●	●

● Single-component pack ▲ Minimum order quantity ■ Can be ordered individually